



# Essential Alpha LED Light

## Technical Specs

Control Unit	Essential Alpha
Lamp weight	11 g
Lamp material	Anodised Aluminium
Operating time	Up to 12 hours
Charging time	Approx. 8 hours
LED diode	Natural white (6,300 Kelvin)
Maximum light intensity at 25 cm	37,000 lux
Light intensity at working distance 35 cm	14,000 lux
Light spot at working distance 35 cm	85 mm
Brightness setting	Stepless
Connection system	Essential Alpha framepart
Cable type	Black super flex
Warranty	Battery and cable 1 year/Lamp 2 years