

Please note that the actual perceived magnification can vary, depending on working distance and prescription.

These examples are to be used as a guide and represent angular magnification, based on a working distance of 45 cm.



No magnification

Galilean HD



Magnification	2.3 x	2.8 x	3.3 x	3.8 x
Field of View	12 cm	10.5 cm	8 cm	6.5 cm
Depth of Field up to	28 cm	20 cm	15 cm	10 cm

Essential



Magnification	2.5 x	3.0 x
Field of View	11 cm	8.3 cm
Depth of Field up to	25 cm	15 cm