



Total Pure & Total Intense LED Light

Technical Specifications

Material	Aluminium
Compatible Control Unit	PowerGo
Connection type	USB-C
Light spot at working distance 35 cm	70 mm uniform light spot
Lamp weight with framepart	18 g
Light colour	Total Pure: 5,000 K Total Intense: 6,500 K
Maximum light intensity at 25 cm	Total Pure: Min. 21,900 Lux Total Intense: Min. 51,500 Lux
Maximum light intensity at 35 cm	Total Pure: Min. 11,300 Lux Total Intense: Min. 25,400 Lux
Colour rendering index	Total Pure: >90CRI Total Intense: >70CRI
LED chip	Total Pure: Samsung LH351D 5000K Total Intense: Osram GW CSSRM3.PM 6500K
Warranty	3 years